



VALCOR SCIENTIFIC
 2 LAWRENCE ROAD
 SPRINGFIELD, NJ 07081
 USA

SV690C Series Operator and Maintenance Manual

S:\TriContinent\VSF-E-0038.docx

REVISIONS

ECO #	REV LTR	PAGE NUMBERS REVISED	DATE	APPROVED
6218	-	Initial Release	11/12/13	

(REVISIONS ARE DETAILED ON THE FOLLOWING PAGE.) NOTICE: Unauthorized use, manufacture or reproduction in whole or in part is prohibited. Drawing, design and other disclosures are property of Valcor Engineering Corporation®. Copyright © 2013, Valcor Engineering Corporation®.

	NAME	DATE
WRITTEN BY:	P. F. Worst	11/12/13
CHECKED BY:		
ENGINEER:		
Q.C.:		
PRODUCTION:		
ENG. MANAGER		

CODE IDENT #	DOCUMENT NUMBER
96487	VSF-E-0038 SV690C OMM SHEET 1 (0580-0040)

TABLE OF CONTENTS

<u>Para.</u>	<u>Title</u>	<u>Sheet</u>
1.	INTRODUCTION.....	4
2.	HARDWARE.....	8
3.	SOFTWARE.....	14
4.	MAINTENANCE.....	23
5.	ORDERING INFORMATION.....	26


1. INTRODUCTION – SV690C SERIES PUMPS

The SV690C Series is a precision OEM pump module designed for use with automated instrumentation. All typical aspiration and dispensing operations expected from an OEM syringe module can be accomplished with this unit. The stepper motor driven syringe will accurately and precisely handle fluids from 5 μ L to 12.5 mL. Standard resolution is 3,000 steps per full syringe stroke and can be set into Hi-Res mode providing 24,000 steps.

A SV690C Series is ideal for experienced instrument manufacturers who are looking to maximize the investment they have already made in existing instrumentation, by using the SV690C Series as a direct replacement for the competitive product. This pump utilizes standard communication protocols and mounting configurations used by other compact pump designs, thereby eliminating the need to perform expensive and time consuming re-development.

This manual was created specifically for the use of instrument designers with the knowledge of existing firmware of similar designs or those looking to upgrade and /or retrofit existing applications. In order to make your product transition quick and simple, a QUICK START CONFIGURATION GUIDE is included in this manual that illustrates only those areas that are unique to the Valcor SV690C Series.

For those customers that are just being introduced to the SV690C Series design, this manual includes basic information you will require to get started. However, a complete SV690C Series Precision Pump Software Manual is available upon request.



SV690C Precision Pump Module

QUICK START CONFIGURATION GUIDE (for replacing an XP3000)

Configuration Jumpers

The SV690C Series has the configuration jumpers located in different locations (see figure on Sheet 7). The functions of the jumpers are also slightly different. Please see below for jumper settings compared to XP3000 with standard software and microstep-enabled software:

Function	SV690C Series	XP (std)	XP (micro)	Action
Protocol ¹	N/A	JP1-2	N/A	Installed = OEM
Reserved ²	J2-1			
EEPROM Mode	J2-2	JP1-3	JP1-3	Installed = Autostart
Overload Detection	J2-3	JP1-1	JP1-1	Installed = Disabled
Communication Rate	J2-4	JP1-4	JP1-4	Installed = 38.4 K/125K
Valve ³	J2-5	JP4, 1-2	N/A	Installed = 120° Valves

Note 1 – The SV690C Series has automatic protocol detection.

Note 2 – A spare jumper is provided on J2-1.

Note 3 – The SV690C Series has a jumper to select between 120° 3-port valve, and all other valve types (90°).

RS485 Termination Jumpers

The SV690C Series has the configuration jumpers located in different locations (see figure Sheet 7). Pumps are shipped with termination jumpers installed. Remove them if not needed.

Function	SV690C Series	XP3000	Action
RS485A	J9-1	JP2, 1-2	Installed = Terminated
RS485B	J9-2	JP2, 3-4	Installed = Terminated
CAN	J9-3	N/A	Installed = Terminated

Software Compatibility


The SV690C Series is 100% compatible with the XP3000 protocols and documented command set. There are some commands and responses that are different as shown below:

Command	Value	Description	Difference
^	< n > 0..255	Sets threshold value for fluid detection	N/A for SV690C ¹
z	None	Set position count to encoder count	N/A for SV690C ¹
?4		Report actual position of plunger	N/A for SV690C ¹
?22		Report current value from fluid sensor	N/A for SV690C ¹
&		Report firmware version	New response ²
#		Report firmware checksum	New response ²

Note 1 – The SV690C Series does not include a valve fluid detector circuit or a high resolution encoder. These commands are accepted, but are not functional.

Note 2 – The SV690C series has different electronics and firmware. These report commands will respond with the same amount of characters, but with different response data.

QUICK START CONFIGURATION GUIDE SV690C SERIES JUMPER LOCATIONS



SPECIFICATIONS

Downloadable specifications can be found on the Valcor website: www.valcor.com

HARDWARE

J3 Mating Connectors

Manufacturer	Description	Manufacturer's P/N
AMP	15 pin female-solder cup, receptacle	747909-2
Cinch	15 pin female-solder cup, receptacle	DA-155
AMP	Plastic shield with male screw retainers	207908-4
Cinch	Plastic shield with male screw retainers	SDH-15GL-CS

J3 Wiring

Power and communication is supplied via a single cable to the SV690C Series.

Pin	Function	Details
1	24 VDC	See Below
2	RS232 TxD line	Output Data
3	RS232 RxD line	Input Data
4	RTS	
5	CAN HI Signal line	
6	CAN LO Signal line	
7	AUX Input #1	TTL Level (4.7K Ω pullup)
8	AUX Input #2	TTL Level (4.7K Ω pullup)
9	Ground	Power and Logic
10	Ground	Power and Logic
11	RS485 A line	
12	RS485 B line	
13	AUX Output #1	TTL Level (1K Ω pullup Max Source/Sink current = ± 20 mA)
14	AUX Output #2	TTL Level (1K Ω pullup Max Source/Sink current = ± 20 mA)
15	AUX Output #3	TTL Level (1K Ω pullup Max Source/Sink current = ± 20 mA)

CAUTION: Do not plug or unplug DB-15 connector with power applied or damage may result.

Power Supply Requirements


Linear or switching supplies may be used.

Output Voltage	24 VDC Nominal $\pm 10\%$
Output Current	≥ 1.5 Amp Peak

Fluid Connections

The SV690C Series has 1/4"-28 fittings as the standard configuration. M-6 connections can also be obtained through special order. Please contact Valcor for details.

ELECTRICAL CONNECTION LOCATIONS SV690C SERIES



CONFIGURATIN JUMPER AND SWITCH SETTINGS

J2

Jumper	Function	Jumper Installed	Jumper Removed
J2-1	Unused	Spare	
J2-2	EEProm Autostart	Enabled	Disabled
J2-3	Plunger Overload	Disabled	Enabled
J2-4	Communication Rate	38,400 baud (RS232/RS485) 125 kbits/s (CAN)	9,600 baud (RS232/RS485) 100 kbits/s (CAN)
J2-5	Valve	120° 3-ports	Other Valves (90° ports)

J9

Jumper	Function	Jumper Installed	Jumper Removed
J9-1	RS485 Termination A	Terminated	Not Terminated
J9-2	RS485 Termination B	Terminated	Not Terminated
J9-3	CAN Termination	Terminated	Not Terminated


Address Switch

The address switch is used to provide each SV690C Series a unique address in a multi-pump application. The address switch has 16 locations (0 – E). The F address location is used to activate the self-test. The address locations can be set by using a small flathead screwdriver and rotating the switch in either direction to the appropriate address.




Address Switch

RS232 CABLING DIAGRAM




* See Sheet 10 for location of J9

RS485 CABLING DIAGRAM



See Sheet 10 for location of J9

CAN CABLING DIAGRAM



SOFTWARE RS/232 AND RS485 COMMAND SUMMARY

Control Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
R	N/A	N/A		Execute command string	
x	0...3	N/A		Execute next command in buffer based on auxiliary inputs	<0> = Execute if Input 2 is Low and Input 1 is Low <1> = Execute if Input 2 is Low and Input 1 is High <2> = Execute if Input 2 is High and Input 1 is Low <3> = Execute if Input 2 is High and Input 1 is High
X	N/A	N/A		Re-execute last executed command string	
G<n>	0...30000	0		Repeat command sequence	0 = Loop forever
g	N/A	N/A		Mark the start of a repeat sequence	
M<n>	0...30000			Delay command execution	Milliseconds
H<n>	0 - 2	0		Halt command execution	<0> = Wait for [R] or either input 1 or 2 to go low <1> = Wait for [R] or input 1 to go low <2> = Wait for [R] or input 2 to go low
T	N/A	N/A		Terminate command	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Initialization Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
Z<n>	0...40	0		Initialize plunger, valve to the right	<0> = Initialize at full plunger force <1> = Initialize at low plunger force <2> = Initialize at lower plunger force than <3> = Initialize at speed 16, low plunger force <4> = Initialize at speed 18, low plunger force <5-9> = same as <0> <10-40> = initialize at defined plunger speed
Y<n>	0...40	0		Initialize plunger, valve to the left	Same as Z<n>
W<n>	0...40	0		Initialize plunger without valve	Same as Z<n>
z<n>	0...3000 (0...24000)	0		Set pump's internal position counter to value specified	
k<n>	0...120 (0...960)		24	Syringe dead volume command	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Plunger Movement Commands (Note: Argument limits are for 3000 step/stroke pumps)

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
A<n>	0...3000 (0...24000)	0		Move plunger to absolute position	
a<n>	0...3000 (0...24000)	0		Move plunger to absolute position, not busy	
P<n>	0...3000 (0...24000)	0		Relative pickup	
p<n>	0...3000 (0...24000)	0		Relative pickup, not busy	
D<n>	0...3000 (0...24000)	0		Relative dispense	
d<n>	0...3000 (0...24000)	0		Relative dispense, not busy	

Valve Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
I	N/A	N/A		Move valve to Input position	
O	N/A	N/A		Move valve to Output position	
B	N/A	N/A		Move valve to Bypass position	
E	N/A	N/A		Move valve to Extra position. Valid for 4 position valves only.	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Set Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
K<n>	0...100 (0...800)		10	Backlash steps	
L<n>	1...20		14	Set acceleration/deceleration slope	
v<n>	0...1000		900	Set start velocity in Hz	
V<n>	1...6000		1400	Set top velocity in Hz	
S<n>	0...40		11	Set speed	
c<n>	1...2700		900	Set cutoff velocity in Hz	
h<n>	0...100		10	Set syringe hold current in %	
m<n>	0...100		75	Set syringe run current in %	
C<n>	0...025		0	Cutoff velocity in steps	
N<n>	0...2		0	Set microstep positioning and velocity mode	<0> = Both microstep position and velocity mode off <1> = Microstep position mode on, velocity mode off <2> = Both microstep position and velocity mode on
J<n>	0...7			Sets the 3 TTL auxiliary outputs	
J<ppppn>	<pppp> 1-3000 <n> 0...7			Set auxiliary outputs based on syringe position	
^<n>	0...255			Non-functional command	
b	N/A	N/A		Non-functional command	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Eeprom Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
s<n>	0...14			Load program string into EEprom	
e<n>	1...14			Execute EEprom string	
U				Set pump configuration parameters	Refer to Software Manual
u	1...16			Set system configuration parameter into EEprom	For factory use only

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Report Commands

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
Q	N/A	N/A		Report system status	
?	N/A	N/A		Report absolute plunger position	
?0	N/A	N/A		Same as ?	
?1	N/A	N/A		Report start velocity in Hz	
?2	N/A	N/A		Report peak velocity in Hz	
?3	N/A	N/A		Report cutoff velocity in Hz	
?4	N/A	N/A		Same as ?	
?5	N/A	N/A		Same as ?	
?6	N/A	N/A		Reports valve position (I, O, B, E)	
?10	N/A	N/A		Report command buffer status	
?12	N/A	N/A		Report number of backlash steps	
?13	N/A	N/A		Report status of Aux 1 input	
?14	N/A	N/A		Report status of Aux 2 input	
?15	N/A	N/A		Report number of pump initializations. Counter not implemented, will always report 1	
?16	N/A	N/A		Report number of plunger movements. Counter not implemented, will always report 1	
?17	N/A	N/A		Report number of valve movements. Counter not implemented, will always report 1	
?18	N/A	N/A		Number of valve movements since last ?18	
?19	N/A	N/A		Reports if pump is initialized with a Y, Z, z, or W command. 0 = not initialized 1 = initialized	
?20	N/A	N/A		Report firmware checksum	
?22	N/A	N/A		Non-functional command for compatibility, will always return 255	
?23	N/A	N/A		Report firmware version	
?24	N/A	N/A		Report the syringe's dead volume as set by the k command	
?25	N/A	N/A		Report hold current in %	
?26	N/A	N/A		Report run current in %	
?27	N/A	N/A		Report EEPROM configuration data as set by the u command	
?28	N/A	N/A		Report if 3 or 4 position valve is installed (jumper status)	
?29	N/A	N/A		Same as Q	
?30-?44	N/A	N/A		Report user EEPROM execution strings	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE
RS/232 AND RS485 COMMAND SUMMARY - Continued

Report Commands - continued

Command	Operand Range <n>* ()=Microstep Mode Limits	Default Operand	Power Up Default	Command Description	Operand Description
F	N/A	N/A		Same as ?10	
&	N/A	N/A		Same as ?23	
#	N/A	N/A		Same as ?20	
RZ	N/A	N/A		Same as ?	
RV	N/A	N/A		Same as & and ?23	
%	N/A	N/A		Same as ?18	

* Operand range listed is for 3000 step/stroke pumps only

SOFTWARE CAN COMMAND SUMMARY

On-the-fly Commands Frame Type = 0

Command	Operands	Command Description
V	Same as RS232/RS485	Top Velocity
T	N/A	Terminate

Action Commands Frame Type = 1

Command	Operands	Command Description
		All RS232/RS485 commands, with the exception of Report commands, are valid Action commands in CAN bus Mode.

Common Commands Frame Type = 2

Command	Operands	Command Description
0	N/A	Reset mode
1	N/A	Start loaded command
2	N/A	Clear loaded command
3	N/A	Repeat last command, like X
4	N/A	Stop action immediately, same as T command

SOFTWARE CAN COMMAND SUMMARY - Continued

Report Commands Frame Type = 4

Command	Operands	Command Description
0	N/A	Plunger position
1	N/A	Reports encoder position, like ?4
2	N/A	Same as report command 0
3	N/A	Reports valve position, like ?6
4	N/A	Top velocity, like ?2
6	N/A	Start velocity, like ?1
7	N/A	Cutoff velocity, like ?3
10	N/A	Buffer status, like F
12	N/A	Backlash steps, like ?12
13	N/A	Input 1 status, like ?13
14	N/A	Input 2 status, like ?14
15	N/A	Number of pump initializations, like ?15. Note, currently not implemented, always returns a 1
16	N/A	Number of plunger movements, like ?16. Note, currently not implemented, always returns a 1
17	N/A	Number of valve movements, like?17. Note, currently not implemented, always returns a 1
18	N/A	Number of valve movements since last report, like ?18
19	N/A	Report if pump is initialized. 1 + initialized, 0 = not initialized
20	N/A	Firmware checksum, like ?20
22	N/A	Non-functional command to maintain backward firmware compatibility. Will always return 255
23	N/A	Firmware version, like &
24	N/A	Syringe dead volume, like ?24
29	N/A	Current status, like Q

MAINTENANCE

Recommended Maintenance

Daily Maintenance	Action
Inspect syringe seals and valves for leaks and proper operation	Replace as required
Inspect tubing fittings for leaks	Tighten or replace as required
Inspect for any fluid or material on outside of pump	Clean as required
"Park" inactive syringes	Flush with DI water at end of use and "park" syringe full of system fluid in full down position


SYRINGE REPLACEMENT

To replace the syringe, prime all fluid from the syringe, and follow the steps below:

1. To replace the syringe, lower the plunger driver.
 - For 3000 step/stroke pumps, if the power is on, this can be done by sending the “A3000R” command or if the pump is powered off, manually lower the plunger by pushing firmly down on the syringe plunger holder.
 - For 24000 step/stroke pumps, with the power on, send the “A24000R” command. You cannot manually lower the plunger on 24000 step/stroke pumps.
2. Turn the plunger lock screw counterclockwise 3 turns.
3. Turn the syringe counter clockwise and remove from the valve.

To install the syringe follow the steps below:

4. Screw the threaded portion of the syringe clockwise into the valve.
5. Fully seat the plunger into the syringe plunger holder.
6. Securely tighten (clockwise) the plunger lock screw.
7. Initialize the pump.



VALVE REPLACEMENT


NOTE: Valcor syringe pump valves are interchangeable among pumps. Refer to the “Configuration Jumper and Switch settings”, (Sheet 11) for instructions.

To replace the valve, prime all fluid from the syringe, and follow the steps below:

1. Remove syringe as described in steps 1 – 3 on Sheet 25.
2. Disconnect the tubing from the valve.
3. Remove the two valve screws used to secure the valve to the panel.
4. Remove the valve from the pump.

To install the valve follow the below steps:

5. Align internal valve “D” slot with the flat on the motor shaft.
6. Gently push and rotate the valve into the correct position to align the pins with the holes in the panel.
7. Securely tighten (clockwise) the plunger lock screws.
8. Initialize the pump.



ORDERING INFORMATION

ACCESSORIES AND REPLACEMENTS

To view or download Accessories and/or Replacement items available for the SV690C Series pumps, please visit our website at www.valcor.com.

ORDERING INFORMATION/CUSTOMER SERVICE

Valcor is committed to exceeding our customers expectations when it comes to our products and services. To place an order or inquire on the SV690C Series or any Valcor product, please contact Valcor at one of the following:

Valcor Engineering Corporation
2 Lawrence Road
Springfield, NJ 07081
(973) 467-8400
(973) 467-9592 (FAX)
scientific@valcor.com
www.valcor.com

WARRANTY AND RETURNS

We design and manufacture our syringe pumps to be the most reliable product available. We stand behind the SV690C Series pumps with a 1 year warranty on material and workmanship, excluding abuse and misuse. If notified of such defects during the warranty period, Valcor will, at its option, either repair or replace products which prove to be defective.

Please refer to the Valcor website www.valcor.com for detailed warranty and return information.